

Type	Material
<b>Dental Alloy</b>	
Non-precious Alloy for PFM & Implant Framework & primary telescopes & secondary telescopes (Co-Cr casting)	Wirobond SG,Bego,DE,CE:0197
Non-precious Alloy for PFM & Implant Framework & primary telescopes & secondary telescopes (Co-Cr laser printing)	Argeloy NP Supreme CE:0086
Non-precious Alloy for Partial Framework (Co-Cr)	Wironit,Bego,DE,CE:0197,:
Titanium	Titan Ti 1/22, DeguDent DE, CE:0124
Argident 87	Argident 87( Au 87.42% Pt 4.5% Pd 5.9% Ag 1% In, Sn, Ir, Fe) Argen CE:0086
Argident 75	Argident 75 (75.1% Au, Pd 12%, Ag 10%,In 1.8%,Sn 1%) Argen CE:0086
Argident 52	Argident 52 ( Au 52.5%,Pd 26.9%, Ag16.0%, In 2.5% ) Argen CE:0086
Argistar 43+	Argistar 43+ (Au 2%, Pd 42.9%, Ag 46%) Argen CE:0086
Argident Euro	Argident Euro (Au 40%, Pd 39.4%, Ag 10%) Argen CE:0086
Argenco W	Argenco W (Pd 24.97% Ag 70.0% In 3.0%) Argen CE:0086
Argenco 75	Argenco 75 (Au 74.5%, Pd 3.5%, Ag 11%, Cu 10.45%) Argen CE:0086
Argenco 50	Argenco 50 (Au- 50% Pd 5.0% Ag 32.5% Cu 11.45% ) Argen CE:0086
<b>Zirconia &amp; e.max &amp; PMMA</b>	
E.max Press	IPS e.max Press, Ivoclar Vivadent, CE:0123
Zirconia	ZexxZr T, Sagemax Boceramics,inc,CE:0086
Full Zirconia	ZexxZr T, Sagemax Boceramics,inc,CE:0086
Temporary Crown	PMMA Disk, YAMAHACHI DENTAL MFG., Co., CE:0120
<b>Ceramic Powder</b>	
Ceramic for PFM	IPS InLine, Ivoclar Vivadent, CE:0123
Ceramic for Titanium	Triceram, DENTAURUM DE,CE:0483
Ceramic for e.max	IPS e.max Ceram, Ivoclar Vivadent, CE:0123
Ceramic for Zirconia	Love Ceram,Dentsply, CE:0124
<b>Composite</b>	
Composite	SR Adoro, Ivoclar Vivadent, CE:0123
<b>Removable Denture</b>	
Resin Teeth	NEW ACE (Ant) YAMAHACHI(Japan), CE:0120
Resin Teeth	NAPERCE (Post) YAMAHACHI(Japan), CE:0120
Acrylic Denture	Mega Press NV(Mega Press NV+Jet), Megadental GmbH,CE:0297
Flexible Denture	TCS, Thermoplastic Comfort Systems,Inc,CE:0050
<b>Orthodontics Material</b>	
Night-guard (Soft-hard)	Erkoflex, Erkodent DE, CE: This material is Class 2 material, doesn't need to have CE No.It can be used under EC regulation.
Night-guard (Pressed, Hard)	Erkodur, Erkodent DE, CE: This material is Class 1 material, doesn't need to have CE No.It can be used under EC regulation.
Bleaching-guard	Erkoflex, Erkodent DE, CE: This material is Class 1 material, doesn't need to have CE No.It can be used under EC regulation.
Anti-snoring	Erkoflex, Erkodent DE, CE: This material is Class 1 material, doesn't need to have CE No.It can be used under EC regulation.